OPTI Rapid Lysis Buffer

English Version

Reagent for the enzymatic release of nucleic acid from upper respiratory samples

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OPTI Rapid Lysis Buffer

Name and Intended Use

OPTI Rapid Lysis Buffer is intended for the rapid release of nucleic acid (DNA and RNA) from upper respiratory specimens such as nasopharyngeal, nasal, and oropharyngeal swab specimens collected in Universal Transport Medium, Viral Transport Medium or saline without nucleic acid extraction. The process involves enzymatic digestion of the sample, followed by heat-inactivation of the proteolytic enzyme. The lysed sample may be used in the OPTI SARS-CoV-2 RT-PCR Test which is validated for qualitative detection of RNA from the SARS-CoV-2 virus.

IMPORTANT: Do not use OPTI Rapid Lysis Buffer on respiratory specimens collected in guanidine based transport medium or in devices that contain DNA/RNA stabilizing agents such as lithium salts.

Note: Do not use OPTI Rapid Lysis Buffer on specimens that appear bloody or turbid. Consider using the OPTI DNA/RNA Magnetic Bead Kit (REF 99-58015) instead.

Materials and Storage

Identification/ General Information	Quantity	Storage		Freeze/Thaw cycles
	500 samples	At Receipt	After Reconstitution	
OPTI Rapid Lysis Concentrate, Lyophilized REF 45-56372-00 Reconstitute to 500 μ L and dilute to final volume of 11 mL in OPTI Rapid Lysis Diluent (see instructions below). The expiration date on the vial is valid for either the dry or reconstituted form.	500 μ L	-25 to 8°C	-25 to -15°C	\leq 6
OPTI Rapid Lysis Diluent REF 45-56376-00 Used for the reconstitution and dilution of the OPTI Rapid Lysis Concentrate (see instructions below).	11 mL	-25 to 8°C		N/A

Note: See table at the end of the insert for a description of symbols used on the insert and labels.

Materials Required but Not Provided

- Pipettes (adjustable volume)
- Nuclease-free, aerosol-resistant pipette tips
- Heating blocks or water baths that can maintain a temperature of 60°C (\pm 2°C) and 97°C (\pm 2°C), accommodate reaction tubes and fit inside a Class II Biological Safety Cabinet.
- 2 mL thermostable reaction tubes with safe lock
- Microcentrifuge for 2 mL microtubes capable of reaching 1500 – 3000 x g
- Vortex or equivalent

Laboratory Practices and Warnings

- Do not use reagents past expiration date.
- Wear buttoned lab coats, powder free gloves, face protection when working with specimens. Refer to Interim Laboratory Biosafety Guidelines for Handling and Processing Specimens Associated with Coronavirus Disease 2019 (COVID-19) from CDC website.
- Use appropriate reaction tubes (e.g. thermostable plastic) that can be securely closed during heating.
- To avoid cross-contamination, use nuclease-free, aerosol-resistant pipette tips, and physically separate the workplaces for nucleic acid extraction/handling, PCR setup and PCR.
- Clinical samples should always be regarded as potentially infectious material.

Preparation of OPTI Rapid Lysis Buffer

OPTI Rapid Lysis Buffer is prepared by mixing Rapid Lysis Concentrate (RL Concentrate) with the Rapid Lysis Diluent (RL Diluent).

- Reconstitute the RL Concentrate to 500 μ L with RL Diluent
- Allow the RL Concentrate to sit at 18 to 26°C for 10–20 minutes
- Mix the RL Concentrate and transfer the solution back into the RL Diluent bottle (final volume is 11 mL)
- Mix the RL Diluent bottle by inversion and change the bottle name to Rapid Lysis Buffer
- Aliquot the Rapid Lysis Buffer as appropriate and store solutions frozen (\leq 6 freeze/thaw cycles)

Procedure

-
- 1 Label each reaction tube according to the specimen ID. Include negative controls according to laboratory protocol. Pipet 20 μ L of OPTI RealPCR Rapid Lysis Buffer to each tube.
-
- 2 Inside a Class II Biological Safety Cabinet, pipet 20 μ L nasopharyngeal, nasal, or oropharyngeal specimen into each tube from step 1. Take precaution to avoid cross-contamination and minimize exposure to the aerosol associated with the clinical specimens.
-
- 3 Close the cap of each tube securely. Vortex briefly.
-
- 4 Incubate for 30 minutes \pm 2 minutes at 60 °C \pm 2°C in a heating block or water bath.
-
- 5 Incubate for 10 minutes \pm 1 minute at 97°C \pm 2°C in a heating block or water bath.
-
- 6 Spin the tubes briefly in a microcentrifuge to draw condensation back down into the tube.
-
- 7 Allow the tubes to cool for 5 minutes at 2 to 8°C or on ice before proceeding to real-time PCR testing. If the sample is not tested immediately, see the "Sample Storage" recommendations below.

Sample Storage

Following the lysis protocol, lysates may be stored up to 4 weeks at –25 to –15°C.

Note: if samples are stored frozen, heat at 60°C (1–3 minutes) or thaw at 18 to 26°C; vortex and microcentrifuge briefly.

GHS Classification

Component	GHS Symbol	Hazard phrases	Precaution phrases
RL Concentrate		Danger 334	261, 285, 304+341, 342+311

Hazard phrases

H 334 May cause allergy or asthma symptoms or breathing difficulties if inhaled.

Precaution phrases

P 261 Avoid breathing mist/vapors.

P 285 In case of inadequate ventilation wear respiratory protection.

P 304+341 If inhaled: If breathing is difficult, remove person to fresh air and keep comfortable for breathing.

P 342+311 If experiencing respiratory symptoms: Call a poison center/doctor.

For further information, please see Material Safety Data Sheets.

For technical assistance on OPTI Rapid Lysis Buffer:

IDEXX USA Tel: +1 800 548 9997 or +1 207 556 4895

IDEXX Europe Tel: +800 727 43399

Contact your IDEXX area manager or distributor or visit our website.

IDEXX and OPTI are trademarks or registered trademarks of IDEXX Laboratories or its affiliates in the United States and/or other countries. All other product and company names and logos are trademarks of their respective holders.

©2021 OPTI Medical is an IDEXX Laboratories, Inc. company. All rights reserved. Made in France

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©2021 OPTI Medical ist ein Unternehmen von IDEXX Laboratories, Inc. Alle Rechte vorbehalten.
OPTI Medical Systems, Inc.
MT Promedt Consulting GmbH